


## BALL-HEAD

The Champions® Implant Ball-Head is a screw-shaped one-piece implant that is inserted to bone level. The implant has diameters of 2.5/3.0/4.0 mm. It is available in lengths of 8/10/12/14/16/18/20/22/24 mm. The implant covers a wide range of indications. The implant is suitable for the insertion in the long-term healed jaw bone (delayed implantation), also for the delayed insertion (6-8 weeks after the tooth extraction) and under certain conditions for immediate implantation (implantation directly after the tooth extraction). The Champions® Implant Ball-Head can be placed in bone in the maxilla and mandible (D1-D4) to support a removable prosthesis. There are at least 3 implants/teeth per quadrant in the maxilla and 2 implants/teeth per quadrant in the mandible.

Please note the restrictions of indications in the Instructions for Use.

Material: Titanium grade 4



| Diameter | Length | Height | Art. No. | Type |
|----------|--------|--------|----------|------|
| 2.5 mm   | 10 mm  | 6.3 mm | 1900     | ●    |
|          | 12 mm  | 6.3 mm | 2015     |      |
|          | 14 mm  | 6.3 mm | 2035     |      |
|          | 16 mm  | 6.3 mm | 2055     |      |
| 3.0 mm   | 8 mm   | 6.8 mm | 1890     |      |
|          | 10 mm  | 6.8 mm | 1905     |      |
|          | 12 mm  | 6.8 mm | 2016     |      |
|          | 14 mm  | 6.8 mm | 2040     |      |
|          | 16 mm  | 6.8 mm | 2060     |      |
|          | 18 mm  | 6.8 mm | 2080     |      |
|          | 20 mm  | 6.8 mm | 2100     |      |
|          | 22 mm  | 6.8 mm | 2120     |      |
| 4.0 mm   | 8 mm   | 7.1 mm | 1895     |      |
|          | 10 mm  | 7.1 mm | 2005     |      |
|          | 12 mm  | 7.1 mm | 2025     |      |
|          | 14 mm  | 7.1 mm | 2050     |      |
|          | 16 mm  | 7.1 mm | 2070     |      |
|          | 18 mm  | 7.1 mm | 2090     |      |
|          | 20 mm  | 7.1 mm | 2110     |      |
|          | 22 mm  | 7.1 mm | 2130     |      |
| 24 mm    | 7.1 mm | 2150   |          |      |