


SQUARE "NEW ART"

The Champions® Implant Square is a screw-shaped one-piece implant that is inserted to bone level. The implants are available in diameters of 3.0/3.5/4.5/5.5 mm and lengths of 6/8/10/12/14/16/18/20/22/24 mm. The implants cover a wide range of indications. The implant is suitable for the insertion in the long-term healed jaw bone (delayed implantation), also for the delayed insertion (6-8 weeks after the tooth extraction) and under certain conditions for immediate implantation (implantation directly after the tooth extraction). The Champions® Implant Square can be placed in bone in the maxilla and mandible (D1-D4).

Please note the restrictions of indications in the Instructions for Use.

Material: Titanium grade 4




Diameter	Length	Height	Art. No.	Type
3.0 mm	10 mm	8.0 mm	1015	○
	12 mm	8.0 mm	1040	
	14 mm	8.0 mm	1080	
	16 mm	8.0 mm	1120	
3.5 mm	6 mm	8.1 mm	0900	
	8 mm	8.1 mm	1000	
	10 mm	8.1 mm	1020	
	12 mm	8.1 mm	1050	
	14 mm	8.1 mm	1090	
	16 mm	8.1 mm	1130	
	18 mm	8.1 mm	1160	
	20 mm	8.1 mm	1170	
	22 mm	8.1 mm	1180	
24 mm	8.1 mm	1190		
4.5 mm	6 mm	8.6 mm	0910	
	8 mm	8.6 mm	1010	
	10 mm	8.6 mm	1030	
	12 mm	8.6 mm	1060	
	14 mm	8.6 mm	1100	
	16 mm	8.6 mm	1140	
5.5 mm	6 mm	9.2 mm	0920	
	8 mm	9.2 mm	1012	
	10 mm	9.2 mm	1035	
	12 mm	9.2 mm	1070	
	14 mm	9.2 mm	1110	
	16 mm	9.2 mm	1150	

SQUARE "CLASSIC"

The Champions® Implant Square Classic is a screw-shaped one-piece implant that is inserted to bone level. The implant has a diameter of 3.5 mm. It is available in lengths of 8/10/12/14/16 mm. The implant covers a wide range of indications. The implant is suitable for the insertion in the long-term healed jaw bone (delayed implantation), also for the delayed insertion (6-8 weeks after the tooth extraction) and under certain conditions for immediate implantation (implantation directly after the tooth extraction). The Champions® Implant Square Classic can be placed in bone in the mandible (D1).

Please note the restrictions of indications in the Instructions for Use.

Material: Titanium grade 4



Diameter	Length	Height	Art. No.	Type
3.5 mm	8 mm	8.1 mm	1005	●
	10 mm	8.1 mm	1025	
	12 mm	8.1 mm	1055	
	14 mm	8.1 mm	1095	
	16 mm	8.1 mm	1135	