

C-SAILOR PRO IMPLANT MOTOR

NEW

The implant motor allows for maximum freedom while working due to its light motor in the handpiece, the easy-to-use touch screen, and the built-in water supply.

INTRODUCTION PRICE

incl. green handpiece
with LED-lighting

€ 1,495
+ VAT



Implant Motor

- 1 x Implant Motor with Touch Screen
- 1 x Handpiece (green 20:1) with LED-Lighting and Tray
- 1 x Multi-Functional Pedal

Delivery in a high-quality aluminium box

REF 6117



Technical Parameters:

- Max. power: 150 W
- Power supply: AC 230V 50 / 60 Hz
- Motor speed: 300-4,000 rpm
- Adjustable torque: 5 Ncm-80 Ncm
- Water supply volume: adjustable in 4 blocks, up to 150 ml/min
- Lighting: up to 30,000 Lux
- Motor: brushless surgical motor
- Dimensions (cm): 42 x 38 x 26
- Total weight: 5.5 kg

Highlights:

- Large touch screen for best visibility and selection of treatment parameters
- Easy handling and installing
- Rinsing device with cooling function
- Built-in water supply, adjustable to 4 blocks
- New motor generation with a torque of up to 80 Ncm
- High-performance, light brushless motor for powerful torque output without power loss or fatigue signs for the user
- Sterilization of motor and cable
- LED-lighting up to 30,000 Lux
- Adjustable lighting and sound
- Multi-functional pedal controller

You can find new
products here:

(scan QR-Code)



Order in the Champions-
Implants-Shop or per Fax!

Further information about our products:

www.champions-implants.com

champions  implants

IMPLANT MOTOR

Order

_____ piece(s) C-Sailor Pro Implant Motor (REF 6117)

€ 1,495*

Billing address (eventually with stamp)

Name _____

Street _____

ZIP Code/City/Country _____

Phone _____ Fax _____

Customer No. _____

Shipping address (eventually with stamp)

Name _____

Street _____

ZIP Code/City/Country _____

Phone _____ Fax _____

Customer No. _____

* Prices + VAT

Fax-Order: +49 (0) 6734 10 53

Free shipping, EU-countries only, for orders from € 500 (net price) up

By placing an order for the following references, the customer accepts the General Terms and Conditions (GTC) of Champions-Implants GmbH (GTC on the homepage > Shop > GTC).

Please send me the current product catalog

Please send me the course program

CHAMPIONS-IMPLANTS GMBH

Im Baumfeld 30 | Champions Platz 1 | D-55237 Flonheim

fon +49 (0)6734 91 40 80 | fax +49 (0)6734 10 53

info@champions-implants.com

Join our courses offered by Future Dental Academy!

Friday: MIMI® – Minimally Invasive Method of Implantation

Saturday: BioWin! zircon implant & Smart Grinder

Day course (Fri. or Sat.) =10 education points!

You can find dates and registration forms on our website:

www.champions-implants.com



Bestellschein Chirurgie Einheit en
18. April 2019, 8:39 vorm.

BE PART OF THE CHAMPIONS® (R)EVOLUTION!

